



COMBAT LITTER STAND

Unlike many predecessors, the FMS Compact Litter Stand can fold to approximately 1/2 its height. This simple feature allows our litter stand to be packed for mobile operations offering advanced treatment in forward areas.



1 - Folded configuration (stowed)	2 - Pull legs apart to deploy	3 - Deployed configuration (low pos.)	4 - Fully deployed configuration
Height 438 mm	Height 427 mm	Height 424 mm	Height 810 mm
Width 116 mm	Width 330 mm	Width 1139 mm	Width 504 mm
Length 733 mm	Length 733 mm	Length 733 mm	Length 733 mm
Cubed 0.0368 cubic m	Cubed 0.103 cubic m	Cubed 0.354 cubic m	Cubed 0.298 cubic m
Weight 4.4 kg	Weight 4.4 kg	Weight 4.4 kg	Weight 4.4 kg

Extreme temperatures affect litter stand surface.
Use care for patient comfort.



Weight tolerance in fully deployed configuration (4).
Tested in pairs.

PRODUCT INFO

- Stands fold small enough for Load Out bags and Hard Cases
- Load tested in pairs to 500 lbs.
- Accessories increase treatment capabilities
- Fielded by USMC (AMAL)
- NSN 6530-01-618-7419
- Globally Sourced / Trade Compliant
- Lightweight and rugged Provides immediate trauma / triage station anywhere
- Ideal for FOB's or shelters
- The CTS and sustainment parts are now available for Direct Purchase by DoD, US.Govt, Federal and Local Law Enforcement
- Originally designed, manufactured and tested by Ferno Military Systems